



AZ Instrument Corp.

Web site: <http://www.az-instrument.com.tw>
E-mail: info@az-instrument.com.tw

Tel: 886-4-25326668 (rep.)
Fax: 886-4-25326593

Edit by Naima 21 Apr. 2020

Question

Can I use Low Ionic Strength pH pen to measure higher Ionic strength solution

Reply:

Yes, AZ's Low Ionic Strength pH pen can also do general pH measurement in the pH range from 0.0~14.0.

Generally, the pH measurement in low ionic strength solution will be unstable, but the performance of AZ's Low Ionic Strength pH pen will remain the same no matter its low ionic strength or normal ionic strength solution for your reference.