



**AZ Instrument Corp.**

Web site: <http://www.az-instrument.com.tw>  
E-mail: [info@az-instrument.com.tw](mailto:info@az-instrument.com.tw)

Tel: 886-4-25326668 (rep.)  
Fax: 886-4-25326593

---

Edit by Naima 9 Apr. 2020

## **Question**

**Does skin color affect readings of 8877?**

## **Reply:**

For model #8877, the infrared emissivity from dark-skinned to light-skinned remains essentially constant (0.98) and does not affect readings for your further reference.