



Edit by Naima 11 May 2020

Question

What is the fever temperature by 8877?

Reply:

Generally, forehead thermometer is a quick screening tool to check if there's any people who's potentially in fever.

To confirm if one is in fever or not, it still need to be tested by doctor.

As far as we know, in most people, a fever is an oral temperature above 38°C. Medical research hasn't determined an exact correlation between oral, rectal, ear, armpit, and forehead temperature measurements, nevertheless, base on the correlation of temperature results, forehead temperature is usually around 0.3°C to 0.6°C lower than an oral temperature.

Here you can find more detail information:

<https://www.mottchildren.org/health-library/tw9223>